

HaloTherapyPRO

Setup & Operating Instructions

Standard Controller

Rev. October 2026

SYSTEM SETUP

System Setup

Step 1 - Unpack the Box

Carefully remove all items from the packaging.

You should have: • HaloTherapyPRO unit • Salt cover (salt lid) • Fan housing with fan installed • Power cord • 1 lb bag of HaloTherapy Salt

WARNING: If anything is missing or damaged, please contact HaloTherapyPRO customer support before use.

Step 2 - Position the HaloTherapyPRO

Place the unit in its designated location.

- Place the HaloTherapyPRO in its designated location, such as a booth, room, or tent.
- Set it on a stable, flat surface with adequate airflow around the unit.

Step 3 - Install the Fan Housing

Connect and seat the fan housing on top of the unit.

- Plug the fan into its connector on the HaloTherapyPRO, if it is not already connected.
- Insert the fan housing into the slot on top of the HaloTherapyPRO.

TIP: Make sure the fan label is facing forward — this indicates the correct direction of airflow.

Step 4 - Add the HaloTherapy Salt

Load the salt bowl with the included HaloTherapy Salt.

- Ensure the salt bowl is properly seated in the unit.
- Place the salt lid over the salt bowl.
- Add 1 teaspoon of the included HaloTherapy Salt.

Step 5 - Power On

Connect power and switch the unit on.

- Plug the included power cord into the HaloTherapyPRO.
- Plug the other end into a suitable wall outlet.
- Switch the unit on.

TIP: That's it — your HaloTherapyPRO is ready to use.

RUN TIME & PAUSE SETTINGS

Operating Cycle Information

Each device is tested and pre-programmed before shipping.

Your HaloTherapyPRO operates on the following pre-set cycle:

- 10 seconds ON
- 30 seconds OFF

Depending on the size of your booth, room, or tent, runtime may be adjusted at your discretion for the ideal session.

TIP: You can adjust the on/off cycle timing to suit the size of your space.

Recommended Settings

- Grind: 10 seconds
- Pause: 30 seconds

TIP: The Grind setting shown may be too high for smaller setups such as pop-up tents or compact salt spaces. Smaller enclosed areas can reach higher salt concentrations more quickly, so you may need to reduce the Grind time (RED) and/or increase the pause time (BLUE) to decrease salt output and find the setting that works best for your setup. Adjustments may vary depending on room size and preference.

Adjusting Time

Adjusting Grind & Pause Times

Change how long your HaloTherapyPRO Grinds and Pauses.

1. Press SET once. The red display (Grind Time) will flash.
2. Use the ▲ / ▼ buttons to set the desired run (grind) time.
3. Press SET again. The blue display (Pause Time) will flash.
4. Use the ▲ / ▼ buttons to set the desired pause time.
5. Wait 6 seconds. The settings will automatically save.

Video tutorial available at HaloTherapyPRO.com/setup

MODE SETUP

Before You Begin — How the Controls Work

Learn how the controls work before starting Step 1.

Before You Begin

Watch the video above to learn how the settings work. Do not change anything yet. When you are ready to program the control panel, begin with Step 1.

P0 and P1 are setting names and do not change.

Hold SET for 3 seconds to enter the settings.

Press SET to switch between P0 and P1.

Use ▲ / ▼ to change the number shown below P0 or P1.

Ready to program? Start with Step 1 below.

Video tutorial available at HaloTherapyPRO.com/setup

Halogenerator Control Panel Setup Guide

Your HaloTherapyPRO arrives **fully programmed and ready to use.** Only follow these instructions if your control panel settings need to be reset.

Control Panel Overview

- * **P0 / 0** = Timer set to seconds
- * **P1 / 5** = Automatic cycle setting
- * **Red display** = Grinder ON time
- * **Blue display** = Grinder PAUSE time

Important: P0 and P1 are the **setting names**. The number shown below each setting is the value that you change.

Step 1 - Set Timer Mode (P0)

Set the bottom value to 0. The P0 setting name will stay the same.

1. **Press and hold SET for 3 seconds** to enter Mode Setup.
2. The display will show **P0**. **P0** is the setting name and will not change.
3. Use the **▲ / ▼** buttons to change the **bottom value** until it shows **0**.
4. The display should now show:
P0
0
5. Press **SET** to continue.



Step 2 - Set Operating Mode (P1)

Set the bottom value to 5. The P1 setting name will stay the same.

****Set the bottom value to 5. The P1 setting name will stay the same.****

****IMPORTANT:**** Complete this step within ****3 seconds**** of Step 1. If you wait too long, the settings will save automatically and you will need to start again at Step 1.

1. After completing Step 1, ****P1**** will appear on the display.
2. ****P1** is the setting name and will not change.**
3. Use the ****▲ / ▼** buttons** to change the ****bottom value**** until it shows ****5****.
4. The display should now show:
****P1****
****5****
5. Wait a few seconds for the settings to ****save automatically****.



TIP: Wait 6 seconds to save settings.

Support & Warranty

For additional information, customer support, or warranty claims, please visit www.HaloTherapyPRO.com. Thank you for choosing HaloTherapyPRO. Enjoy your halotherapy experience.